

LOW LEVEL LIGHT
CARBON INFRARED HEAT
COMBINATION THERAPY

Prophylaxis

Therapy

Regeneration

Training


Sport

Competition



EN

FOR HORSES



MOYAVE

H45

Tendon boot ensures perfect fit

Carbon infrared heat
with 8 to 12 μ m

45 LEDs with 632 nm



Inner surface

Breathable functional material

Operated via Control Box **CB22**
Long-life integrated battery






FOR HORSES



FOR DOGS

MOYAVE

P45



All-purpose flat pad



45 LEDs with 632 nm

Carbon infrared heat
with 8 to 12 μ m

2 Velcro® fasteners



60 cm

Suitable for many
parts of the body



Operated via Control Box **CB22**
Long-life integrated battery





Your veterinarian or therapist can advise you about a possible therapy and treatment

MUSCLES AND FASCIA

Promotes circulation
Reduces cramps
Muscle disorders
Pulled/strained muscles and muscle fiber rupture

TENDONS AND JOINTS

Bone spavin
Excessive strain
Insertion desmopathies
Periosteal inflammation
General swelling and inflammation
Inflammation of a tendon and its sheaths


NERVES

Localized nerve damage

PAIN RELIEF

Spondylosis, Dysplasia of the hip
Kissing Spine Syndrome

**POSSIBLE
APPLICATIONS**



PRODUCTS

Delivery includes:



Control Box **CB22**



Carrying Case **CC1**



Mains Adapter **PA12**

MOYAVE **H45**
TENDON BOOT LEFT



MOYAVE **H45**
TENDON BOOT RIGHT



MOYAVE **P45**
FLAT PAD



2 Velcro® fasteners
60 cm

All these products are available from our online store: www.moyave.com





TECHNOLOGY and BENEFITS

A COMBINATION OF TECHNOLOGIES

Low-level light therapy with red light at 632 nm (LED)

Long-wave infrared heat at 8 to 12 μm (FIR)

Multi-layer technology with breathable functional materials



LOW-LEVEL-LIGHT
45 HIGH POWER LED's



CARBON INFRARED HEAT
240 cm² TREATMENT AREA

BENEFITS

Safe and straightforward handling


Robust and high-quality finish

Flexible application for all sizes

Self-contained thanks to Li-Ion battery operation

Nine preset treatment programs

LED and FIR can also be used on their own



A PRACTICAL
EXAMPLE

Erhard Schulze - veterinarian specializing in acupuncture and gold acupuncture with horses and small animals

Thanks to the MOYAVE Pad, I at last have a device I can use to make therapeutic and prophylactic interventions. It's effective, sits nicely in my hand, and is easy to clean. The combination of red LED light and infrared heat means we can achieve good results. I have a small animal practice, and I like to use the pad after administering gold acupuncture to dogs. We've had some excel-

lent experiences in our practice, particularly in the case of spondylosis patients where I like to place the pad on their back during the waking-up phase. The pad really relaxes the musculature, and this means that the overall procedure causes very little distress to the animal. The preset programs make it easy to cater precisely to the dog's individual needs.



Monika Bartholomé - veterinarian specializing in acupuncture for horses

Our horses often suffer from blockages in the back, the lumbar region, or the sacroiliac joint. This is where the pad really comes into its own. Handling is very straightforward, and the horses are very willing to accept it.

We can also use the MOYAVE pads and tendon boots in all sorts of ways for indications such as scar disruption, swollen tendons, or cases of tendinitis.



MOYAVE® - a FROBAS brand
FROBAS GmbH - Gebrüder-Eicher-Ring 45 - 85659 Forstern - Germany
Phone: +49 8124 91890 40 - info@moyave.com - www.moyave.com



www.facebook.com/moyave